



Transparent Budget Cryoboxes in Polypropylene - resistant down to -196C



Dimensions: 50 x 130 x 130 mm and 80 x 130 x 130 mm
Holds 81 vials, 9 x 9
Suited for 1.2ml and 2ml cryogenic vials.
Usable temperature range of -196°C to +121°C
Printed lid with codes for easy identification.
Five colours; red, yellow, green, clear and blue

Cryoboxes, 50mm height

Cat. No.	Description	Pack size	Colour	Price in sFr.
02-5000-8100	81-Place Cryobox for 1,2 and 2ml tubes, with clear lid, 50mm high	10	Natural	59.00
02-5000-8101	81-Place Cryobox for 1,2 and 2ml tubes, with clear lid, 50mm high	10	Blue	59.00
02-5000-8102	81-Place Cryobox for 1,2 and 2ml tubes, with clear lid, 50mm high	10	Green	59.00
02-5000-8104	81-Place Cryobox for 1,2 and 2ml tubes, with clear lid, 50mm high	10	Red	59.00
02-5000-8106	81-Place Cryobox for 1,2 and 2ml tubes, with clear lid, 50mm high	10	Yellow	59.00
02-5000-8199	81-Place Cryobox for 1,2 and 2ml tubes, with clear lid, 50mm high	10	Mixed colours	59.00

Prices exclude MWSt

Freight and packing costs to be charged for small orders under sFr. 300.-; all other orders free delivery.

Life Systems Design AG

Chilefeldstrasse 17
 CH-5634 Merenschwand

Tel. 056/664 59 80, Fax. 056/664 59 88

E-Mail: mail@lifesystemsdesign.ch

Web-site: www.lifesystemsdesign.ch